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First read the whole instructions and adjust the surface of the paper wick accordingly.

## 1. Construction

300ml formic acid-bottle with filling, emptying scale and childproof top

Dropper

Perforated paper-wick

Baseplate with fixing prongs

## 2. Mode of operation

The paper-wick is placed on the baseplate. The bottle filled with FORMIVAR® is inverted and placed on the baseplate. The formic acid slowly drops through the dropper of the bottle and is absorbed by the paper-wick. From the wet paper-wick, formic acid evaporates within the hive. The formic acid concentration of the hive air increases slowly but steadily. As soon as the paper-wick is fully saturated the formic acid concentration in the hive remains constant. Thus the Liebig-Dispenser allows an even and bee-friendly formic acid treatment from which the bees experience no problems.

## 3. Instruction of use

### 3.1 Safe use of formic acid

Formic acid is highly corrosive. It has to be handled with great care. In any case **protective glasses, acid-proof gloves and long-sleeved clothes should be worn.** Splashes on the skin have to be washed off thoroughly at once. Have water ready! If in contact with the eyes wash thoroughly with water. Consult a doctor immediately!

### 3.2 Filling of formic acid into the bottle

The bottles are filled at a suitable place with formic acid. For this the dropper has to be removed from the neck of the bottle (e.g. with small flat pliers). The necessary amount of formic acid (see table 1) is measured with the help of the filling scale on the bottle. The dropper is then

replaced and the bottle is secured with its childproof top for transport.

## 3.3 Adjustment of the paper-wick area for the regulation of evaporation rate

### 3.3.1 Base adjustment

The ability to change the **size** of the paper-wick makes it possible to regulate the required rate of evaporation to cater for a multitude of influences (temperature, kind and **size** of hive, outside temperature, sun strength). The base adjustment to the area of the paper-wick should serve as a guide (see table last page). The paper-wick area can simply be reduced by tearing along the perforations.

### 3.3.2 Adjustment during the treatment

It is strongly recommended that you note the rate of evaporation in the initial days of the treatment and keep to the guidelines for the rate of evaporation according to **table 1** (last column).

- Immediately after insertion of the Liebig-Dispenser the acid level in the bottle sinks relatively quickly. Thus the initial level should only be read and noted after 1 hour for the control of the evaporation rate.
- It is essential to include weather prediction in the decision.

#### The evaporation rate is too high:

Reduce the paper-wick area by tearing along the perforations. **Attention:** It is essential that safety measures (chapter 3.1) are observed when the wick is already soaked with formic acid.

#### The evaporation rate is too low:

Enlarge the paper-wick area again by replacing the torn-off pieces. Replaced pieces should overlap slightly with the paper-wick that is already soaked with formic acid.

## 3.4 Carrying out the treatment

**Caution:** The colonies should stand as level as possible. Otherwise there is the danger that the formic acid could drip from the baseplate into the colony.

At high temperature or when using cold formic acid, an overpressure may develop in the bottle causing a leakage

**Table 1:** dosage of formic acid for 1<sup>st</sup> and 2<sup>nd</sup> treatment as well as daily evaporation rate in ml.

Kind of hive	1 <sup>st</sup> treatment (July/August)		2 <sup>nd</sup> treatment (September)		Evaporation guidelines per colony and day	
	Amount of FORMIVAR®		Amount of FORMIVAR®			
	85%*	60%*	85%*	60%*	85%	60%
One chamber hive	50ml	90 ml	100ml	140 ml	10–15 ml	20–25 ml
Two chamber hive	100ml	140 ml	200ml	200 ml	20–30ml	30–40 ml
Dadant, DNM	100ml	140 ml	200ml	200 ml	15–20 ml	25–30 ml
Warm way	50ml	90 ml	100ml	140 ml	10–15 ml	20–25 ml

\*For technical reasons, the imprint of the graduation of the bottle may vary +/- 10 ml. However this small inaccuracy has no influence on the success of the treatment.

of formic acid into the bee colony. Therefore, **only formic acid of ambient temperature should be used**. Treatment is preferably started in the early morning hours or during a period with moderate temperatures but **not at hot summer conditions**.

1. For getting enough place for the Liebig-Dispenser, add an empty super.
2. Shut any open floors. Check if hive and lid are closed well.
3. Reduce paper-wick if necessary by tearing off pieces (see last page).
4. Put adapted paper-wick over fixing prongs on the baseplate. (Attention: Use only one paper-wick.)
5. Place the baseplate with the paper-wicks in the centre of the brood frames. The paper-wick should cover thereby as much as possible of the occupied combs.
6. Put the formic acid bottle with dropper head first onto the fixing prongs of the baseplate. For improved stability of the bottle on the baseplate, a piece of paper (approx. 3 cm of kitchen paper or paper-wick) may be added between bottle and fixing prongs.
7. Close the colony. Keep the entrance of the hive normally open.
8. Monitor the daily evaporation rate (guidelines see table 1, last column).
9. Remove the Liebig-Dispenser as soon as the formic acid bottle is empty (after one to two weeks).

### 3.5 Cleaning

Clean the base plate at room temperature. Remove propolis with rubbing alcohol or methylated spirits. Remove residual wax with an oven-foam cleaner or with warm water (max. 40° C) and washing-up liquid.

## 4. Concept of control

### 4.1 Time

The Liebig-Dispenser is used exclusively after the honey-flow in late summer. Two treatments are sufficient:

- **First treatment** lasting 3–6 days; **before** feeding (middle July–start August).
- **Second treatment** lasting 6–14 days; **after** feeding (September).

### 4.2 Winter treatment


Winter treatment with oxalic acid in the brood-less time may only be omitted, if the natural mite-fall in the middle of november is **under 1.0 mite per day on average**.

If, in the coming year, drone brood is cut out or nuclei are made, only those colonies with more than 1 naturally fallen mites per day have to have the follow-up treatment. The methods for using oxalic acid are described on page 11.

### 5. Possible side effects

The Liebig-Dispenser has proved itself as bee-friendly in both one and two chamber Zander hives as well as Warm way. Young brood and mature brood are especially sensitive to formic acid, therefore an effective formic acid treatment always leads to a certain loss of brood. These losses have however, as found by research, no negative influence on the consequent development of the colony.

**Queen losses** can be excluded if the **guidelines for rate of evaporation are not exceeded**.

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## English

### Treatment with oxalic acid

Oxalic acid treatments guarantee a constant and high efficacy against the varroa mite in broodless colonies (November to December). These treatments are also well tolerated by the bees.

#### Vaporising oxalic acid

At outdoor temperatures of at least 4°C, 1–2 g of oxalic acid dihydrate is vaporised through the entrance of the hive with the VARROX®-vaporiser.

#### Safety precautions:

Rubber gloves, Protective mask (type EN 149: 2001 FFP3).

#### Trickling oxalic acid

At outdoor temperatures of at least 3°C, 5 ml of lukewarm solution of oxalic acid sugar water (about 25°C) is trickled directly on the bees in every occupied beeway.

#### Quantity per hive:

- 30 ml weak colony
- 40 ml normal colony
- 50 ml strong colony

#### Safety precautions:

Protective glasses, rubber gloves.






#### Spraying oxalic acid

At outdoor temperatures of at least 8°C, 3–4 ml of a 2.1% solution of oxalic acid is sprayed per side of comb occupied with bees.

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